

## Elite™ Homogeneous Mitochondrial Membrane Potential Assay Kit

CATALOG NUMBER: CA-M185, 10 plates

### Description

Mitochondrial membrane potential is an important parameter of mitochondrial function used as an indicator of cell health. Elite™ Homogeneous Mitochondrial Membrane Potential Assay Kit uses a fluorescent cationic, lipophilic dye (**m-MPI**) to signal the loss of mitochondrial potential. This kit is the non-wash version. In normal or healthy cells, **m-MPI** concentrates in the mitochondrial matrix where it forms red fluorescent aggregates. However, in apoptotic and necrotic cells, **m-MPI** diffuses out of mitochondria and changes to green fluorescent monomeric form. This kit provides all the necessary reagents for the analysis of mitochondrial integrity in cells. Fluorescence excitation/emission maxima: 514/529 nm for monomer form; 585/590 nm for J-aggregate form.

### Features

- **Solubility:** Excellent water solubility, no precipitation unlike widely used JC-1 dye.
- **Specificity:** Selectively entering into mitochondria.
- **Convenient:** non-wash
- **Application Platform:** Fluorescent microplate reader, flow cytometry, and fluorescence microscopy.
- **Robust and Highly Sensitive:** Excellent response to membrane potential change (Green/Red ratio). m-MPI signal changes proportionally to mitochondrial permeability defined by electrochemical gradient ( $\Delta\Psi$ ) across the membrane.

### Applications

- Cell apoptosis study.
- Cell toxicity.
- High Throughput Screening of apoptosis activator or inhibitors.

### Kit Components

- **Component A:** Mitochondrial Membrane Potential Indicator (m-MPI) in DMSO 200  $\mu$ L (1000x)
- **Component B:** m-MPI Assay Buffer 15 mL (10x)
- **Component C:** Signal enhancer 1 mL (200x)

### Storage

Mitochondrial Membrane Potential Indicator (m-MPI) should be protected from light and stored at -20 °C.

**Note: Aliquot Component A into 10 vials (20  $\mu$ l each) and store at -20 °C when received.**

Signal enhancer should be stored at room temp.

m-MPI Assay Buffer should be stored at 4°C.

### Materials Required (but not supplied)

- 96- or 384-well microplates: Tissue culture microplate with black wall and clear bottom is recommended
- Fluorescence microplate reader or flow cytometer

### ASSAY PROTOCOL

#### 1. Cell Preparation

1.1 Harvest cells when they reach 80-90% confluence in flasks. Trypsinize cells, and count a portion of the cells with a hemocytometer.

**Note: It is very important that the cells DO NOT reach >90% confluence.**



- 1.2 Add 20  $\mu$ L/well of cell suspension to 384-well plates (black/clear) the day before the experiment. The cell number needs to be optimized for each assay. For assays performed on fluorescence plate reader, optimal assay conditions require a confluent monolayer of cells prior to the assay (6K cells/well with HeLa cells and 8K cells/well with HepG2 cells).
- 1.3 Allow cells to attach by leaving the cell plates at room temperature for 30 minutes. Transfer the plates to a cell culture incubator and grow the cells overnight.

## 2. Preparation of 2X Dye-Loading Solution

- 2.1 Remove a vial of aliquoted m-MPI (**Component A**) from  $-20^{\circ}\text{C}$ , and allow it to thaw at room temperature.
- 2.2 Prepare 1X m-MPI assay buffer by adding 135 mL ddH<sub>2</sub>O into 15 mL 10X m-MPI assay buffer (**Component B**).
- 2.3 Preparing 2X Dye-Loading Solution for one cell plate by transferring 20  $\mu$ L of m-MPI Stock Solution (**Component A**) and 100  $\mu$ L signal enhancer (**Component C**) to 10 mL of 1x m-MPI Assay Buffer (from step 2.2). Mix well by vortexing for 10 seconds.
- 2.4 Aliquot the unused m-MPI, and store at  $-20^{\circ}\text{C}$ .

**Note: Do not re-use 2x Dye-loading Solution. Avoid repetitive freeze-thaw of m-MPI DMSO stock.**

## 3. Dye Loading and Mitochondrial Membrane Potential Assay

### 3.1 Kinetic assay with on-line compound addition

- a) Dilute testing compounds in 1x HBSS buffer at 5x final concentrations.
- b) Remove cell plates from incubator; Add 20  $\mu$ L of 2X Dye-Loading Solution to each well and incubate at  $37^{\circ}\text{C}$  for 30 min.
- c) Place the cell plates on a fluorescence plate reader with on-line liquid handling system (e.g. FlexStation or FDSS), and perform the mitochondrial membrane potential assay with on-line addition of 10  $\mu$ L of test compounds to the cell plates.
- d) For assays performed on a FlexStation (MDS), use the following wavelength parameters. Green channel (Monomer form): Excitation: 485 nm; Emission: 530 nm; Cutoff 515 nm; Red channel (J-aggregates): Excitation: 485 nm; Emission: 590 nm; Cutoff 570 nm.

**Note: Dispense speed and height for compound additions need to be optimized for each instrument.**

### 3.2 Endpoint assay

- a) Prepared 5X testing compounds in 1x HBSS.
- b) Remove cell plates from incubator; add 5  $\mu$ L of 5x compound solution. Incubate at  $37^{\circ}\text{C}$  for proper amount of time (For FCCP and oligomycin, 60 min incubation is recommended).
- c) Add 25  $\mu$ L 2X dye-loading solution to each well and incubate at  $37^{\circ}\text{C}$  for 30 min.
- d) Read on a fluorescence plate reader with bottom read-head. Test filter settings to optimize fluorescence light collection and eliminate bleedthrough of excitation light to emission filters.

### 3.3 Assay on BD Pathway 855 (Cat No. 341036) and BD Pathway 435 (Cat No. 641250) Bioimagers

- a) Dye-loading procedure is the same as described above.
- b) Image the plates on a BD Pathway™ Bioimager using appropriate filters (table 1).


**Table 1. Pathway 855 Bioimager Filter Sets**

Channels	Excitation	Emission	Dichroic
Green (Monomer)	488/10	Fura/Fitc	Fura/Fitc
Red (J-aggregates)	548/20	Fura/Tritc	570LP

**Mitochondrial Membrane Potential Quantitation:** The ratio of green (monomer) to red (J-aggregates) fluorescence is used to quantify mitochondrial membrane potential changes.

Data Analysis

Dose Response Curves of FCCP in HepG2 and Hela Cells



**Figure 1.** End point assays performed on FlexStation (MDS). **(A)** HepG2 Cells treated with FCCP for 60 min. **(B)** Human Hepatocytes treated with FCCP for 60 min. Fluorescence signals were recorded by a FlexStation. Dose response curves were plotted by Prism (GraphPad).